

# FASTER ORTHODONTIC TREATMENT

**AcceleDent® Aura** is a simple-to-use, hands-free device designed to accelerate tooth movement **with only 20 minutes daily use.**

## Fast

- » AcceleDent Aura advances the science of orthodontic care with innovative SoftPulse Technology®.
- » Safe acceleration of the bone remodeling process complements conventional orthodontic treatment.

### AcceleDent Aura delivers compelling clinical results

- » Clinically proven to accelerate tooth movement up to 50%.<sup>1</sup>
- » SoftPulse Technology is based on demonstrated scientific methods in orthopedic and orthodontic bone remodeling.<sup>2,3</sup>

## Safe


### Demonstrated safe and reliable in U.S. clinical trials

- » Investigators reported no risk of increased root resorption and no risk in increased loss of posterior anchorage.<sup>1,3</sup>
- » AcceleDent Aura is a Class II medical device cleared by the FDA.

## Gentle

### AcceleDent may improve the orthodontic experience for patients.

- » SoftPulse Technology micropulses are gentle and exert up to 8x less force than a power toothbrush.



*SoftPulse Technology® generates precisely calibrated micropulses to accelerate remodeling activity during orthodontic treatment.*

**AcceleDent®**  
aura

“AcceleDent gives our patients an attractive option to have the smile they desire in a time frame they will accept.”

*Robert Rudman, DDS, MS*

## YOUR ACCELEDENT AURA COMES WITH FOUR ECONOMICALLY FRIENDLY COMPONENTS:



### ACTIVATOR

The heart of AcceleDent Aura, the small, lightweight Activator generates gentle micropulses. It includes a USB interface which can be plugged directly into a computer to view your usage history via the FastTrac Report.



### MOUTHPIECE

Chosen specifically for you by your orthodontist, the Mouthpiece provides a comfortable fit and snaps easily on and off the Activator for transport and cleaning.



### USB EXTENSION CABLE with POWER ADAPTER

Allows for easy, convenient charging.



### TRAVEL CASE

The Travel Case keeps the components safe, clean and secure.

#### References:

1. Data on file.
2. Rubin C, Judex S, Qin YX. Low-level mechanical signals and their potential as a non-pharmacological intervention for osteoporosis. *Age Ageing*. 2006 Sep; 35 suppl 2:ii32-ii36.
3. Kau CH. A radiographic analysis of tooth morphology following the use of a novel cyclical force device in orthodontics. *Head Face Med*. 2011 Aug 9; 7-14.

AcceleDent<sup>®</sup> is an orthodontic accessory intended for use during orthodontic treatment. Refer to the Directions for Use with AcceleDent<sup>®</sup> Aura or [acceleddent.com](http://acceleddent.com) for warnings and precautions prior to using this product. Federal (U.S.) law restricts this device to sale by or on the order of a licensed doctor. Rx use only for a single patient. Made in the USA ©2013 OrthoAccel<sup>®</sup> Technologies, Inc. AcceleDent<sup>®</sup> is a registered trademark of OrthoAccel<sup>®</sup> Technologies, Inc.